

Product Datasheet

Equine Growth Hormone ELISA Kit (orb439752)

Description	Equine Growth Hormone (GH) competitive inhibition ELISA Kit
Reactivity	Equine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01245
Sensitivity	0.085ng/mL
Entrez	100034180
Expiration Date	Please enquire.