

Product Datasheet

Rat Growth Hormone ELISA Kit (orb439757)

| | |
|--------------------------|------------------------------------------------------------------------------------------|
| Description | Rat Growth Hormone (GH) Sandwich ELISA Kit |
| Reactivity | Rat |
| Range | 0.156-10ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma and other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | P01244 |
| Sensitivity | 0.054ng/mL |
| Entrez | 24391 |
| Expiration Date | Please enquire. |