

Product Datasheet

Bovine GHRH ELISA Kit (orb439758)

Description	Bovine Growth Hormone Releasing Hormone (GHRH) competitive inhibition ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P63292
Sensitivity	5.9pg/mL
Entrez	281191
Expiration Date	Please enquire.