

Product Datasheet

Goat GHRH ELISA Kit (orb439759)

Description	Goat Growth Hormone Releasing Hormone (GHRH) Sandwich ELISA Kit
Reactivity	Goat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P63293
Sensitivity	6.4pg/mL
Entrez	102181610
Expiration Date	Please enquire.