

Product Datasheet

Sheep GHRL ELISA Kit (orb439770)

Description	Sheep Ghrelin (GHRL) competitive inhibition ELISA Kit
Reactivity	Sheep
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q30DT1
Sensitivity	0.045ng/mL
Entrez	442996
Expiration Date	Please enquire.