

## Product Datasheet

### Human GLP1 ELISA Kit (orb439795)

<b>Description</b>	Human Glucagon Like Peptide 1 (GLP1) competitive inhibition ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Sensitivity</b>	5.44pg/mL
<b>Expiration Date</b>	Please enquire.