

Product Datasheet

Human GLP1R ELISA Kit (orb439799)

Description	Human Glucagon Like Peptide 1 Receptor (GLP1R) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P43220
Sensitivity	0.118ng/mL
Entrez	2740
Expiration Date	Please enquire.