

Product Datasheet

Mouse GnRH ELISA Kit (orb439827)

Description	Mouse Gonadotropin Releasing Hormone (GnRH) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P13562
Sensitivity	4.33pg/mL
Entrez	14714
Expiration Date	Please enquire.