

Product Datasheet

Mouse GPER ELISA Kit (orb439857)

Description	Mouse G Protein Coupled Estrogen Receptor 1 (GPER) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8BMP4
Sensitivity	0.054ng/mL
Entrez	76854
Expiration Date	Please enquire.