

## Product Datasheet

### Rat GRIA1 ELISA Kit (orb439887)

<b>Description</b>	Rat Glutamate Receptor, Ionotropic, AMPA 1 (GRIA1) Sandwich ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	15.625-1000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<b>P19490</b>
<b>Sensitivity</b>	6.1pg/mL
<b>Entrez</b>	<b>50592</b>
<b>Expiration Date</b>	Please enquire.