

## Product Datasheet

### CD133 antibody (FITC) (orb43989)

|                        |  |
|------------------------|--|
| <b>Description</b>     | CD133 antibody (FITC)  |
| <b>Species/Host</b>    | Rabbit   |
| <b>Reactivity</b>      | Human  |
| <b>Conjugation</b>     | FITC   |
| <b>Immunogen</b>       | Recombinant Human Prominin-1 protein (508-792AA)                 |
| <b>Target</b>          | PROM1  |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b> | Liquid   |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze. |
| <b>Note</b>            | For research use only  |
| <b>Isotype</b>         | IgG  |
| <b>Clonality</b>       | Polyclonal   |
| <b>Uniprot ID</b>      | <a href="#">O43490</a>   |
| <b>Dilution Range</b>  | ELISA  |
| <b>Expiration Date</b> | 12 months from date of receipt.                                  |