

## Product Datasheet

### Human GSTa1 ELISA Kit (orb439924)

<b>Description</b>	Human Glutathione S Transferase Alpha 1 (GSTa1) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625-40ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>P08263</b>
<b>Sensitivity</b>	0.265ng/mL
<b>Entrez</b>	<b>2938</b>
<b>Expiration Date</b>	Please enquire.