

Product Datasheet

Mouse GSTa1 ELISA Kit (orb439925)

Description	Mouse Glutathione S Transferase Alpha 1 (GSTa1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P13745
Sensitivity	0.129ng/mL
Entrez	14857
Expiration Date	Please enquire.