

Product Datasheet

Human GSTm1 ELISA Kit (orb439932)

Description	Human Glutathione S Transferase Mu 1 (GSTm1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P09488
Sensitivity	0.52ng/mL
Entrez	2944
Expiration Date	Please enquire.