

Product Datasheet

Porcine GSTp ELISA Kit (orb439936)

Description	Porcine Glutathione S Transferase Pi (GSTp) Sandwich ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P80031
Sensitivity	22.8pg/mL
Entrez	100739163
Expiration Date	Please enquire.