

Product Datasheet

Human HBg1 ELISA Kit (orb439989)

Description	Human Hemoglobin Gamma 1 (HBg1) competitive inhibition ELISA Kit
Reactivity	Human
Range	468.75-30000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and erythrocyte lysates.
Sample Types	serum, plasma and erythrocyte lysates.
Uniprot ID	P69891
Sensitivity	122.6ng/mL
Entrez	3047
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)