

Product Datasheet

Porcine HIF1a ELISA Kit (orb440038)

| | |
|--------------------------|--|
| Description | Porcine Hypoxia Inducible Factor 1 Alpha (HIF1a) Sandwich ELISA Kit |
| Reactivity | Porcine |
| Range | 15.625-1000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | A0T1H6 |
| Sensitivity | 5.8pg/mL |
| Entrez | 396696 |
| Expiration Date | Please enquire. |