

Product Datasheet

Human HMGA1 ELISA Kit (orb440056)

Description	Human High Mobility Group AT Hook Protein 1 (HMGA1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P17096
Sensitivity	0.124ng/mL
Entrez	3159
Expiration Date	Please enquire.