

Product Datasheet

Mouse Hpt ELISA Kit (orb440084)

Description	Mouse Haptoglobin (Hpt) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q61646
Sensitivity	13.9ng/mL
Entrez	15439
Expiration Date	Please enquire.