

Product Datasheet

Porcine Hpt ELISA Kit (orb440085)

Description	Porcine Haptoglobin (Hpt) Sandwich ELISA Kit
Reactivity	Porcine
Range	78.125-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8SPS7
Sensitivity	25.5ng/mL
Entrez	397061
Expiration Date	Please enquire.