



Product Datasheet

Rad51 antibody (orb44013)


Description	Rad51 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human DNA repair protein RAD51 homolog 1 protein (1-228AA)
Target	RAD51
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q06609
Dilution Range	WB:1:2000-5000, IHC-P:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(lane 1), J...



Immunohistochemical staining of human li...



Immunofluorescence analysis of HeLa cell...