

Product Datasheet

Porcine HSPA1A ELISA Kit (orb440138)

Description	Porcine Heat Shock 70kDa Protein 1A (HSPA1A) Sandwich ELISA Kit
Reactivity	Porcine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6S4N2
Sensitivity	1.14ng/mL
Entrez	396648
Expiration Date	Please enquire.