

Product Datasheet

Human HSPA8 ELISA Kit (orb440146)

Description	Human Heat Shock 70kDa Protein 8 (HSPA8) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P11142
Sensitivity	0.134ng/mL
Entrez	3312
Expiration Date	Please enquire.