

Product Datasheet

Human HTR1A ELISA Kit (orb440175)

Description	Human 5-Hydroxytryptamine Receptor 1A (HTR1A) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P08908
Sensitivity	0.059ng/mL
Entrez	3350
Expiration Date	Please enquire.