

Product Datasheet

Mouse IAPP ELISA Kit (orb440192)

Description	Mouse Islet Amyloid Polypeptide (IAPP) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P12968
Sensitivity	24pg/mL
Entrez	15874
Expiration Date	Please enquire.