

Product Datasheet

Guinea pig INS ELISA Kit (orb440257)

Description	Guinea pig Insulin (INS) competitive inhibition ELISA Kit
Reactivity	Guinea pig
Range	140.62-9000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P01329
Sensitivity	64.3pg/mL
Entrez	100379579
Expiration Date	Please enquire.