

Product Datasheet

Human INS ELISA Kit (orb440258)

Description	Human Insulin (INS) competitive inhibition ELISA Kit
Reactivity	Human
Range	78-5000pg/ml
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, cell culture supernatants and other biological fluids.
Uniprot ID	P01308
Sensitivity	31.7pg/mL
Entrez	3630
Expiration Date	Please enquire.