



Product Datasheet

CD98 antibody (Biotin) (orb44026)

Description	Mouse Monoclonal to CD98.
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, IHC-P, IP
Immunogen	RAJ1 human Burkitt's lymphoma cell line
Target	CD98
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 3 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	Q01650
Entrez	8140
Dilution Range	Flow cytometry: Recommended dilution: 3 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-filled)



Flow cytometry surface staining pattern ...