

Product Datasheet

Rat INSL3 ELISA Kit (orb440261)

Description	Rat Insulin Like Protein 3 (INSL3) competitive inhibition ELISA Kit
Reactivity	Rat
Range	0.234-15ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9WUK0
Sensitivity	0.062ng/mL
Entrez	114215
Expiration Date	Please enquire.