

Product Datasheet

Mouse INSL6 ELISA Kit (orb440264)

Description	Mouse Insulin Like Protein 6 (INSL6) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9QY05
Sensitivity	30pg/mL
Entrez	27356
Expiration Date	Please enquire.