

Product Datasheet

Human IREM1 ELISA Kit (orb440280)

Description	Human Immune Receptor Expressed On Myeloid Cells 1 (IREM1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8TDQ1
Sensitivity	0.062ng/mL
Entrez	146722
Expiration Date	Please enquire.