

Product Datasheet

Mouse ITGa3 ELISA Kit (orb440297)

Description	Mouse Integrin Alpha 3 (ITGa3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q62470
Sensitivity	0.118ng/mL
Entrez	16400
Expiration Date	Please enquire.