

Product Datasheet

Human ITGb3 ELISA Kit (orb440307)

Description	Human Integrin Beta 3 (ITGb3) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05106
Sensitivity	13.3pg/mL
Entrez	3690
Expiration Date	Please enquire.