

Product Datasheet

Mouse ITLN1 ELISA Kit (orb440321)

Description	Mouse Intelectin 1 (ITLN1) Sandwich ELISA Kit
Reactivity	Mouse
Range	2.5-160pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D9J349
Sensitivity	0.9pg/mL
Entrez	16429
Expiration Date	Please enquire.