

Product Datasheet

Rat KISS1 ELISA Kit (orb440368)

Description	Rat Kisspeptin 1 (KISS1) competitive inhibition ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q7TSB7
Sensitivity	8.47pg/mL
Entrez	289023
Expiration Date	Please enquire.