

Product Datasheet

Mouse KISS1R ELISA Kit (orb440370)

Description	Mouse Kisspeptin Receptor (KISS1R) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q91V45
Sensitivity	0.061ng/mL
Entrez	114229
Expiration Date	Please enquire.