

Product Datasheet

Human KLK12 ELISA Kit (orb440383)

Description	Human Kallikrein 12 (KLK12) competitive inhibition ELISA Kit
Reactivity	Human
Range	15.625-1000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9UKRO
Sensitivity	7.14ng/mL
Entrez	43849
Expiration Date	Please enquire.