

Product Datasheet

Mouse KLK1 ELISA Kit (orb440388)

Description	Mouse Kallikrein 1 (KLK1) Sandwich ELISA Kit
Reactivity	Mouse
Range	18.75-1200pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P15947
Sensitivity	7.52pg/mL
Entrez	16612
Expiration Date	Please enquire.