

Product Datasheet

Human KLK5 ELISA Kit (orb440397)

Description	Human Kallikrein 5 (KLK5) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, saliva, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, saliva, tissue homogenates and other biological fluids.
Uniprot ID	Q9Y337
Sensitivity	29pg/mL
Entrez	25818
Expiration Date	Please enquire.