

## Product Datasheet

### Human KRT33A ELISA Kit (orb440432)

<b>Description</b>	Human Keratin 33A (KRT33A) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.125-5000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<b>076009</b>
<b>Sensitivity</b>	29pg/mL
<b>Entrez</b>	<b>3883</b>
<b>Expiration Date</b>	Please enquire.