

Product Datasheet

Human KRT81 ELISA Kit (orb440438)

Description	Human Keratin 81 (KRT81) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q14533
Sensitivity	0.053ng/mL
Entrez	3887
Expiration Date	Please enquire.