

Product Datasheet

Human LACRT ELISA Kit (orb440446)

Description	Human Lacritin (LACRT) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	saliva, tears and other biological fluids.
Sample Types	saliva, tears and other biological fluids.
Uniprot ID	Q9GZZ8
Sensitivity	0.055ng/mL
Entrez	90070
Expiration Date	Please enquire.