

Product Datasheet

Human LCN1 ELISA Kit (orb440478)

Description	Human Lipocalin 1 (LCN1) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, saliva, tears and other biological fluids.
Sample Types	serum, plasma, saliva, tears and other biological fluids.
Uniprot ID	P31025
Sensitivity	29pg/mL
Entrez	3933
Expiration Date	Please enquire.