

Product Datasheet

Canine LEP ELISA Kit (orb440510)

Description	Canine Leptin (LEP) Sandwich ELISA Kit
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	002720
Sensitivity	0.127ng/mL
Entrez	403616
Expiration Date	Please enquire.