

Product Datasheet

Canine LIF ELISA Kit (orb440537)

| | |
|--------------------------|---|
| Description | Canine Leukemia Inhibitory Factor (LIF) Sandwich ELISA Kit |
| Reactivity | Canine |
| Range | 31.25-2000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | J9NRL6 |
| Sensitivity | 14.7pg/mL |
| Entrez | 403449 |
| Expiration Date | Please enquire. |