

Product Datasheet

Mouse LMNA ELISA Kit (orb440564)

Description	Mouse Lamin A/C (LMNA) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P48678
Sensitivity	0.115ng/mL
Entrez	16905
Expiration Date	Please enquire.