

Product Datasheet

Human LRP1 ELISA Kit (orb440617)

Description	Human Low Density Lipoprotein Receptor Related Protein 1 (LRP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q07954
Sensitivity	0.061ng/mL
Entrez	4035
Expiration Date	Please enquire.