

Product Datasheet

Human LRP6 ELISA Kit (orb440622)

Description	Human Low Density Lipoprotein Receptor Related Protein 6 (LRP6) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O75581
Sensitivity	0.059ng/mL
Entrez	4040
Expiration Date	Please enquire.