

Product Datasheet

Human LTA4H ELISA Kit (orb440631)

Description	Human Leukotriene A4 Hydrolase (LTA4H) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P09960
Sensitivity	0.105ng/mL
Entrez	4048
Expiration Date	Please enquire.