

Product Datasheet

Mouse LTF ELISA Kit (orb440636)

Description	Mouse Lactoferrin (LTF) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, breast milk and other biological fluids.
Sample Types	serum, plasma, breast milk and other biological fluids.
Uniprot ID	P08071
Sensitivity	0.241ng/mL
Entrez	17002
Expiration Date	Please enquire.